**FUJIFILM**

**medwork**

## CAN-Series

CAN1-A3

CAN1-A5

CAN1-A6

CAN1-B3

CAN1-B4

CAN1-D5

404450305893

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
## Instructions for use

### Intended use


The CAN series instruments are used for probing the papilla and for filling the bile duct and the pancreatic duct with contrast agent.

### Product features

CAN1-A6-18-200-35

CAN1-B3-18-260-35

### Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

### General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

**ATTENTION!** Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

### Indications

The indications correspond to the general indications for endoscopic retrograde cholangiopancreatography, cholangioscopy and duodenoscopy.

### Contraindications

The contraindications for the CAN series are the same as the specific contraindications for endoscopic retrograde cholangiopancreatography.

### Possible complications

Pancreatitis, cholangitis, aspiration, perforation, haemorrhage, infection, sepsis, allergic

reaction to contrast medium or other medication, hypertension, respiratory depression or arrest, cardiac arrhythmia or arrest.

### Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**Attention!** When using ERCP catheters with stylet cap check that the Luer Lock is correctly connected.

### Instructions for use

Introduce the ERCP-catheter into the endoscope accessory channel in short increments, until the device is endoscopically visible. Following probing, contrast medium may be injected to fluoroscopically confirm the desired position of the device, while the stylus is in place.

The guide wire connection must be closed for application of contrast agent before using the CAN1-A6-18-200-35.

The X-ray mark helps to precisely determine the position of the catheter tip. For the use of dual-lumen ERCP catheters of the CAN series (e.g., CAN1-B4-18-200-35), we recommend the administration of a dilute contrast agent solution. After removing the stylus from the ERCP-catheter, a guide wire may be inserted to protect access. Flush the catheter lumen with sterile saline solution or distilled water before each introduction of the guide wire. A catheter filled with contrast medium remains may make it difficult to advance the guide wire. To facilitate advancement, immerse every guide wire in sterile saline solution prior to each introduction.

### After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

### List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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