PAP-Series
AXS_tome /AXS_tome⁺

PAP1-A1	PAP1-K2
PAP1-A2	PAP1-K8
PAP1-C2	PAP1-M2
PAP1-C5	PAP1-N2
PAP1-C7	PAP1-P8
PAP1-C8	PAP1-Q7
PAP1-D2	PAP1-S8
PAP1-D7	PAP1-T8
PAP1-D8	PAP1-X7
PAP1-F8	PAP1-X8
PAP1-K1	

404450305718

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
Intended use

The PAP series instruments are used for transendoscopic exploration of the papilla, visualisation of the duct systems and papillotomy, and may also be used for injection of contrast media.

Product features

Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

⌚ ATTENTION! Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

Indications

The indications correspond to the general indications for endoscopic retrograde cholangiopancreatography, cholangioscopy and duodenoscopy.

Contraindications

The contraindications for the PAP series are the same as the specific contraindications for endoscopic retrograde cholangiopancreatography.

Possible complications

Pancreatitis, cholangitis, perforation, haemorrhage, aspiration, infection, sepsis, allergic reaction, hypertension, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, tissue damage, obstruction of the pancreatic duct, fulguration and burns, hyperamylasemia, recurrent bile duct stone, papillary stenosis, biliary stricture, cholecystitis.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Attention! AXS_tome+ papillotomes have a stop function that activates at a tension of approx. 90° to avoid excess tension on the cutting wire (> 90°) in use. The 90°-tension of the papillotome can be overcome with a small amount of force by pulling the double ring handle towards the thumb ring. This prevents a tension of > 90° without impairing the function of the papillotome.

Precautions if using high-frequency current

The instruments are suitable for the use of monopolar HF current. Before using the HF control unit, familiarise yourself with the device manufacturer's safety procedures and instructions for use, and in particular check that the earthing and closed circuit (return electrode) are correct. Observe the relevant HF device manufacturer's recommended setting parameters before use. The instruments must be used with a Type BF or CF generator. The peak voltage measured in the "Cutting/Coagulation Mode" is 1,600 V for PAP Series instruments. If this value is exceeded, an accidental current flow may occur. Do not use the instruments near flammable liquids, explosive gases or in oxygen-enriched environments. During HF application, make sure the instruments do not have any contact with metal objects such as clips etc. because the contact may cause functional failure or fracture of instruments.

Instructions for use

Remove the distal shaping wire from the guide wire lumen of the multi-lumen PAP series papillotome. Carefully inspect the connecting cable for kinks or breaks. Connect the return electrode and the papillotome to the HF control unit when the latter is switched off. Make sure the connectors are joined properly. Insert the papillotome into the working channel of the endoscope and push it forward a little at a time until the tip is visible on the endoscopic image. The multi-lumen papillotome can be inserted using a previously positioned guide wire. Flush the lumen with sterile saline solution before insertion. You can check that the tip of the papillotome is in the correct position by application of contrast agent through the Luer lock connector. Check the settings of the HF control unit and switch it on. Perform the endoscopic papillotomy (EPT). You can change the deflection of the distal instrument tip by operation of the dual-ring handle. After the EPT switch off the HF control unit and remove the connection cable from the handle. Slowly pull the instrument out of the endoscope to prevent contamination of the patient and user with blood and other body fluids, and an associated risk of infection.

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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