FUJIFILM

medwork

POL-Series _ CrossSnare

POL1-X1

404450305781

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
Instructions for use

Intended use

The POL Series instruments are used for tissue resection in the upper and lower gastrointestinal tract.

Product features

Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

ATTENTION! Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

Indications

Polyps, unclear gastrointestinal mucosa lesions.

Contraindications

The contraindications for the POL series are the same as the specific contraindications for oesophagogastroduodenoscopy, colonoscopy, sigmoidoscopy and rectoscopy.

Possible complications

Fulguration of the mucous membrane, perforation, haemorrhage, aspiration, infection, fever, sepsis, allergic reaction to contrast agents or medicines, hypertension, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, transmural burns, thermal injuries, explosion, postpolypectomy syndrome, serositis.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, DO NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Attention! Please ensure that the CrossShare is suitable for both tissue resection with a high-frequency (HF) current and mechanical resection. Studies have shown that there is a reasonable risk of bleeding when polyps with a diameter up to 9 mm are resected without HF current/mechanically. Polyps larger than 9 mm should therefore be resected with HF current.

Precautions if using high-frequency current

The instruments are suitable for the use of monopolar HF current. Before using the HF control unit, familiarise yourself with the device manufacturer's safety procedures and instructions for use, and in particular check that the earthing and closed circuit (return electrode) are correct. Observe the relevant HF device manufacturer's recommended setting parameters before use. The instruments must be used with a Type BF or CF generator. The peak voltage measured in the "Cutting/Coagulation Mode" is 1,600 V for POL Series instruments. If this value is exceeded, an accidental current flow may occur. Do not use the instruments near flammable liquids, explosive gases or in oxygen-enriched environments. During HF application, make sure the instruments do not have any contact with metal objects such as clips etc. because the contact may cause functional failure or fracture of instruments.

Instructions for use

Insert the instrument (with a loop in the tube) into the working channel of the endoscope and push it forwards a little at a time until the catheter tip is visible in the endoscopic image. Adjust the endoscope tip on the tissue to be resected and carefully push the loop out of the tube. Place the snare round the tissue to be resected and pull it slowly into the tube by operating the dual-ring handle until you feel resistance.

Only for mechanical resection/ without HF current:

The tissue is resected by being mechanically detached. The tensile strength required for this procedure is applied at the handle by pulling the dual-ring handle towards the thumb ring.

For resection with HF current:

Before switching the HF control unit, please ensure that the loop only goes around the tissue to be resected, because there is otherwise a risk of a mucous membrane fulguration and possible perforation. Connect the neutral electrode and instrument to the HF control unit while it is switched off. Make sure the plug connections fit correctly. Check the setting of the HF control unit and switch it on. Carry out the planned resection. After resection, switch the HF control unit off and take the connection cable off the handgrip. Slowly pull the instrument out of the endoscope to prevent contamination of the patient and the user by blood and other bodily fluids, and an associated risk of infection.

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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