# POL-Series \_ CrossSnare ZERO

POL1-Z2

404450305911

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
## Intended use

The POL series instruments are used for tissue resection in the upper and lower gastrointestinal tract.

## Product features

## Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

## General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

**Ⓢ ATTENTION!** Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the packaging contains tears or perforations, the seal is not intact or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

## Indications

Polyps, unclear gastrointestinal mucosa lesions.

## Contraindications

The contraindications for the POL series are the same as the specific contraindications for oesophagogastroduodenoscopy, colonoscopy, sigmoidoscopy, rectoscopy and enteroscopy.

## Possible complications

Perforation, haemorrhage, aspiration, infection, fever, sepsis, allergic reaction to contrast agents or medicines, hypertension, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, serositis.

## Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**Attention!** Please note that the POL1-Z series snares are suitable exclusively for tissue resection without HF current, i.e. for mechanical resection. Cold removal with a polypectomy snare is recommended for polyp sizes up to 9 mm.

## Instructions for use

Insert the instrument (with a snare in the tube) into the working channel of the endoscope and push it forward a little at a time until the catheter tip is visible in the endoscopic image. Adjust the endoscope tip on the tissue to be resected and carefully push the snare out of the tube by pressing the control. Place the snare around the tissue to be resected and press the control to pull the snare slowly into the tube until you feel resistance.

The tissue is resected by mechanical removal. Manually apply the tensile strength required for this to the handle by pulling the control quickly towards the thumb ring. Slowly pull the instrument out of the endoscope to prevent contamination of the patient and user with blood and other body fluids, and an associated risk of infection.

## After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

## List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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