FUJIFILM

medwork

PRO-Series _ KENA

PRO1-K1

PRO1-K3

PRO1-P1

PRO1-P2

PRO1-S1

PRO3-C2

PRO3-C3

PRO3-T2

PRO3-T3

404450305810

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
Instructions for use

Intended use

The products of the PRO series are used for endoprosthesis insertion for drainage of the biliary system.

Product features

PRO3-C3

* This product feature only applies to the short-wire-capable PRO1-K3, PRO1-S1-09, PRO3-C3, and PRO3-T3 products.

Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

ATTENTION! Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

Indications

Bile duct obstructions, biliary leaks, biliary strictures, therapy of biliary stones.

Contraindications

The contraindications for the PRO series are the same as the specific contraindications for endoscopic retrograde cholangiopancreatography.

Possible complications

Pancreatitis, cholangitis, perforation, haemorrhage, aspiration, infection, sepsis, allergic reaction, hypertension, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, tissue damage, endoprosthesis migration, fistulas, jaundice, stone formation, biliary colic, inflammation, stricture formation, foreign body in patient, obstruction, occlusion, pain, discomfort.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks,

fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Attention! We recommend implanting FUJIFILM medwork endoprostheses with the aid of the FUJIFILM medwork endoprosthesis introducer sets. This ensures optimal matching of all components. Check whether the flaps of the endoprostheses can be fully unfolded. When using the product, observe the guide wire diameter recommended on the product label. If long-term therapy is required, the endoprosthesis should be assessed for replacement at intervals of 3 months. This endoprosthesis is not intended for use as a permanent implant. Proceed with particular caution when removing endoprostheses in order to avoid injuries.

Instructions for use of guide catheter and pusher

When using short wire, insert the placed guide wire at the distal end of the guide catheter, and let it exit from guide catheter through the lateral guide wire port.

Introduce the guide catheter into the endoscope until the distal tip of the guide catheter is visible through the endoscope. After successful probing, contrast medium can be injected through the guide catheter. To inject the contrast medium, use the Luer lock port on the proximal squeeze connector and attach your syringe there. Before injecting the contrast medium, fix the squeeze connector of the guide catheter in place. Load the guide catheter with a suitable endoprosthesis by removing the proximal squeeze connector. Push the pusher and the endoprosthesis into the endoscope while holding the guide catheter stably in place by means of the Albaran lever. Verify the correct position of the endoprosthesis under fluoroscopy and under endoscopic vision. Carefully withdraw the guide wire and the guide catheter (if you are using one). Make sure you do not displace the pusher, as this keeps the endoprosthesis in position. Next, remove the pusher from the endoscope.

Instructions for use of endoprostheses and (preloaded) placement system

Bring the guide wire into position endoscopically above the stricture, the calculus or the injury. Load the introducer system with a suitable endoprosthesis if necessary. Remove the notice sign if necessary, adjust the protruding length of the guide catheter to suit the requirements, and fix the guide catheter in place with respect to the pusher. 7 Fr endoprostheses are brought into position via the guide wire without a guide catheter. Insert the entire introducer system via the guide wire through the working channel of the endoscope.

If you are using a short wire (only possible with the PRO1-K3, PRO1-S1-09, PRO3-C3, and PRO3-T3 products): Pull the mandrel back behind the guide wire port of the placement system (approx. 2 cm) and fix the mandrel in place on the guide catheter via the cover cap. Introduce the guide wire which is lying in place at the distal end of the guide catheter and exit the system at the lateral guide wire port.

If you are using a long wire: Remove the mandrel from the placement system, introduce the guide wire which is lying in place at the distal end of the introducer system and push the guide wire into the introducer system until it comes out at the cover cap.

Working without wire: push the mandrel right up to the cover cap.

In order to facilitate insertion of the endoprosthesis into the working channel of the endoscope, all FUJIFILM medwork endoprosthesis introducer systems are equipped with an insertion guide. As soon as the distal flap of the endoprosthesis enters the working channel, push the insertion guide over the proximal flap until this has also completely disappeared in the working channel. Then withdraw the insertion guide. Take care to ensure that the insertion guide is not completely inserted into the working channel. With the Albaran lever opened, advance the introducer set in short pushes under radiological control until the radiopaque mark on the guide catheter is positioned above the stenosis in order to ensure adequate guidance of the endoprosthesis.

When using short wire (only possible with the PRO1-K3, PRO1-S1-09, PRO3-C3, and PRO3-T3 products): The prosthesis can be repositioned as long as the guide catheter is placed distally on the guide wire. For this, pull back the insertion system placed above the guide wire together with the pre-loaded prosthesis and the guide wire.

Now unfasten the squeeze connector of the pusher. Using the pusher, advance the endoprosthesis via the guide catheter which is fixed in place. Verify correct positioning of the endoprosthesis under fluoroscopy and under endoscopic vision. Carefully withdraw the guide wire and the guide catheter (if you are using one). Make sure you do not displace the pusher, as this keeps the endoprosthesis in position. Next, remove the pusher from the endoscope. Withdraw the instruments from the endoscope slowly in order to prevent contamination of the patient and the users with blood and other body fluids and to avoid the associated risk of infection.

Instructions for use

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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