SPE-Series

SPE1-B5
SPE1-B6
SPE1-C6

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
Intended use

The SPE1-B6- and SPE1-C6-series instruments are used for application of fluids to the mucous membrane of the gastrointestinal tract.

The SPE1-B5-series instruments are used for rinsing and suction in the bronchial tract.

Product features

Qualification of user

The instruments must only be used with comprehensive knowledge of the technical principles of gastro-intestinal or bronchial endoscopes. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

⌚ ATTENTION! Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

Contraindications

The contraindications for the SPE series are the same as the specific contraindications for oesophagogastroduodenoscopy, colonoscopy, sigmoidoscopy and rectoscopy and bronchoscopy.

Possible complications

Perforation, haemorrhage, infection, sepsis, allergic reaction to contrast agents or medicines, hypertension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Instructions for use

Check that the stylet cap is correctly closed on the Luer lock connector of the spray catheter. Insert the spray catheter into the working channel of the endoscope and advance it a little at a time until it is visible at the exit of the endoscope.

After placement of the spray catheter (SPE1-B6 and SPE1-C6), you can apply colouring agent, contrast agent or medical fluids from the Luer Lock port with the stylet horizontal. The stylet cap must be closed tightly for this.

After placement of the spray catheter SPE1-B5-18-120, you can suction fluids from the bronchial tract or rinse it through the Luer-Lock port with the stylet in place. The stylet cap must be tightly closed for this.

After placement of the spray catheter SPE1-B5-21-120, you can suction fluids from the bronchial tract or rinse it through the Luer Lock port after removing the stylet.

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	SBS		Sterile barrier system

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