**FUJIFILM**

**medwork**

## DIL-Series

DIL1-A1

de en fr es it nl pt sv da el tr  
hr cs et hu lv lt pl sl bg ga mt  
ro sk fi no is

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
**Intended use**

The DIL series instruments are used for the endoscopic dilation of strictures in the gastrointestinal tract.

**Product features****Qualification of user**

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

**General instructions**

Use this instrument only and exclusively for the purposes described in these instructions.

 **ATTENTION!** Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

**Indications**

Strictures in the upper and lower gastrointestinal tract, gastroenteric anastomotic strictures, benign biliary strictures, achalasia, Crohn's strictures, biliary and pancreatic stones.

**Contraindications**

The contraindications for the DIL series correspond to the specific contraindications for oesophagogastroduodenoscopy and endoscopic retrograde cholangiopancreatography.

**Possible complications**

Allergic reaction to contrast agents, bleeding, perforation, sepsis/infection, pancreatitis, cholangitis, foreign body in patient, anxiety, nausea, discomfort, aspiration/inhalation, pain, abscess, fever, fistula, hematoma, laceration, cholecystitis, recurrent bile duct stones, hyperamylasemia.

**Precautions**

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**Attention!** Only use sterile water, saline or contrast agent solutions for inflation purposes. Do not use air or gases for inflation purposes, as these reduce the effectiveness of the balloon. Monitor the working pressure during the procedure. The working pressure may not exceed maximum pressure. Both specifications are provided on the label on your product. Only use the dilation balloon in combination with a dilation instrument equipped with an appropriate pressure display. Bear in mind that the balloon can become damaged during dilation if it comes into contact with sharp-edged or metal objects.

**Attention!** Do not expand the balloon with air or liquid before use as this can damage the endoscope. Before use, remove the protective sheath from the distal end, as well as the metal wire from the guide wire lumen. Connect a dilation instrument to the insufflation hub of the dilation balloon. Apply medical grade silicone oil to the balloon to facilitate forward movement of the balloon in the working channel. Before insufflation, the entire dilation balloon must be pushed fully out of the endoscope, fluoroscopically identified and positioned. If you are using a guide wire measuring approx. 200 cm in length, make sure that it exits the guide wire lumen again at the Guide Wire Port. Loosen the guide wire fixture at the biopsy valve and insert the deflated dilation balloon into the working channel. Then secure the guide wire again in the fixture intended for this purpose. If you are working with a guide wire measuring 400 cm or more in length, the guide wire can exit the guide wire lumen at the Guide Wire Port or at the distal end of the guide wire hub. Note the following when using the Guide Wire Port. The proximal black marking indicates that the Guide Wire Port will soon become visible. The fixation must then be loosened and the dilation balloon removed completely from the working channel. Secure the guide wire. Stop pulling back the dilation balloon immediately if you notice any resistance and then remove the instrument from the patient together with the endoscope. Should pressure suddenly drop or the balloon burst, the procedure must be stopped immediately. Deflate the balloon and remove the instrument from the patient together with the endoscope. Then separate the balloon at the black distal marking. The plastic tubing can now be removed from the working channel.

**Instructions for use**

Insert the guide wire at the distal end into the guide wire lumen. Using short nudging movements, insert the instrument carefully up to the exit portal visible through the endoscope and into the required dilation position. Ensure that the balloon has fully exited the endoscope and verify fluoroscopically that the balloon is correctly positioned within the stricture using the radiopaque markings. Now expand the dilation balloon fully to the inflation pressure shown on the display by adding fluid. The relevant working pressure is provided on the label on your product. To remove the instrument, generate negative pressure using the dilation instrument. The dilatation instrument indicator should be within vacuum range. Ensure that the balloon is fully deflated and that the fluid has been removed completely. Slowly remove the fully deflated balloon from the working channel.

## Instructions for use

### After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

### List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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