FUJIFILM

medwork

POL-Series_SnareInflator

POL1-H3

POL1-H4

de en fr es it nl pt sv da el tr
hr cs et hu lv lt pl sl bg ga mt
ro sk fi no is

404450305884

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
Intended use

The POL series instruments are used for tissue resection in the upper and lower gastrointestinal tract.

Product features**Qualification of user**

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

ATTENTION! Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

Indications

Polyps, unclear gastrointestinal mucosa lesions.

Contraindications

The contraindications for the POL series are the same as the specific contraindications for oesophagogastroduodenoscopy, colonoscopy, sigmoidoscopy and rectoscopy.

Possible complications

Perforation, haemorrhage, aspiration, infection, fever, sepsis, allergic reaction to contrast agents or medicines, hypertension, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, transmural burns, thermal injuries, explosion, postpolypectomy syndrome, serositis.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious

incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

ATTENTION! Test the needle-snare combination to make sure that the needle retracts and extends completely. The instrument should not be used in the ileum and in endoscopes with a lateral telescope (Albaran lever) because deflection may cause instrument failure. When injecting substances, observe the warnings and notes on use provided by the respective manufacturer.

Precautions if using high-frequency current

The instruments are suitable for the use of monopolar HF current. Before using the HF control unit, familiarise yourself with the device manufacturer's safety procedures and instructions for use, and in particular check that the earthing and closed circuit (return electrode) are correct. Observe the relevant HF device manufacturer's recommended setting parameters before use. The instruments must be used with a Type BF or CF generator. The peak voltage measured in the "Cutting/Coagulation Mode" is 1,600 V for POL series instruments. If this value is exceeded, an accidental current flow may occur. Do not use the instruments near flammable liquids, explosive gases or in oxygen-enriched environments.

During HF application, make sure the instruments do not have any contact with metal objects such as clips etc. because the contact may cause functional failure or fracture of instruments.

Instructions for use

Connect the return electrode and the polypectomy handle to the HF control unit when the latter is switched off. Make sure the connectors are joined properly. Connect a syringe filled with injectable solution to the Luer lock connector of the injection needle and flush the air out of the catheter. Insert the instrument into the working channel of the endoscope with retracted polypectomy loop and injection needle. The needle is fully retracted into the catheter as soon as the operating element at the handle has latched noticeably when pulling back. Do not pull the operating element back beyond the latching point. Push the catheter forward in small steps until the catheter tip becomes visible in the endoscopic image. Make sure the needle is completely retracted inside the catheter in order to prevent patient injuries and damage to the endoscope. Align the endoscope tip with the tissue to be resected. Extend the needle by pushing the control and perform the injection. After the injection the needle must be fully retracted into the catheter tip. Carefully extend the snare out of the catheter by operating the dual-ring handle. Place the snare around the tissue to be resected and pull it slowly into the catheter until you feel resistance. Before switching on the HF control unit, make sure the snare only goes round the tissue to be resected because otherwise there is a risk of mucous membrane fulguration and possible perforation. Check the setting of the HF control unit and switch it on. Carry out the planned resection. After resection, switch the HF control unit off and remove the connection cable from the handle. Slowly pull the instrument out of the endoscope to prevent contamination of the patient and user with blood and other body fluids, and an associated risk of infection.

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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