
# PRO-Series \_ KENA

PR01-C2  
PR01-T2

404450305811

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
## Intended use

The PRO series products are used for endoprosthesis insertion for drainage of the biliary system.

## Product features

## Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

## General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

**Ⓢ ATTENTION!** Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

## Indications

Bile duct obstructions, biliary leaks, biliary strictures, therapy of biliary stones.

## Contraindications

The contraindications for the PRO series are the same as the specific contraindications for endoscopic retrograde cholangiopancreatography.

## Possible complications

Pancreatitis, cholangitis, perforation, haemorrhage, aspiration, infection, sepsis, allergic reaction, hypertension, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, tissue damage, endoprosthesis migration, fistulas, jaundice, stone formation, biliary colic, inflammation, stricture formation, foreign body in patient, obstruction, occlusion, pain, discomfort.

## Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the compe-

tent authority of the Member State in which the user and/or patient is established.

**Attention!** We recommend implanting FUJIFILM medwork endoprostheses with the aid of the FUJIFILM medwork endoprosthesis introducer sets. This ensures optimal matching of all components. Check whether the flaps of the endoprostheses can be fully unfolded. If long-term therapy with endoprostheses is required, the endoprosthesis should be assessed for replacement at intervals of 3 months. This endoprosthesis is not intended for use as a permanent implant. Proceed with particular caution when removing endoprostheses in order to avoid injuries.

## Instructions for use

Place the guide wire endoscopically above the stricture, the calculus or the injury. Load the introducer system with a suitable endoprosthesis if applicable. 7 F endoprostheses are brought into position via the guide wire without guide catheter. Introduce the entire introducer system via the guide wire through the working channel of the endoscope. **Attention!** In order to facilitate introduction of the endoprosthesis into the working channel of the endoscope, all FUJIFILM medwork endoprosthesis introducer systems are equipped with an insertion aid. As soon as the front flap of the endoprosthesis enters the working channel, push the insertion aid over the rear flap until this has also completely disappeared in the working channel, then withdraw the insertion aid. Take care to ensure that the insertion aid is not completely inserted into the working channel. With Albarran lever opened, advance the introducer system in short pushes under radiological control, until the radiopaque marking of the guide catheter is positioned as far as possible above the stenosis in order to ensure adequate guidance of the endoprosthesis. Now, disengage the Luer-Lock connection on the guide catheter. Advance the endoprosthesis by means of the pusher via the guide catheter which is fixed in place. Verify the correct position of the endoprosthesis under fluoroscopy and under endoscopic view. Carefully withdraw the guide wire and the guide catheter (if such is present). Make sure you keep the pusher fixed in place, as the pusher keeps the prosthesis in position. Then remove the pusher from the endoscope.

## After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

## List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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