**FUJIFILM**

**medwork**

**clipmaster. -Series**

SPE1-X1-26-220

SPE1-X3-26-220

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
## Instructions for use

### Intended use

The instruments of the SPE1-X-series consist of single-use instruments for the endoscopic application of clips in the gastrointestinal tract.

### Product features

### Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

### General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

**ATTENTION!** Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

### Indications

Mucosal and submucosal defects, ulcer bleeding, Mallory-Weiss syndrome, arteries < 2mm, polyps < 15mm, endoscopical marking.

### Contraindications

The contraindications for the Clipmaster series are the same as the specific contraindications for oesophagogastroduodenoscopy, colonoscopy, sigmoidoscopy, rectoscopy and endoscopic retrograde cholangiopancreatography.

### Possible complications

Infection, haemorrhage, allergic reactions on nickel, swelling, edema, perforation, laceration, foreign body in patient.

### Precautions


In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.


Check the distal end of the spiral. The first clip must still be completely inside the clamp, and it must not extend partly or completely out of the clamp. Only use the Clipmaster if the haemostasis can be checked visually on the endoscopic image after application. The Clipmaster should not be used in endoscopes with side-view optics because the bending of the spiral above the Albaran lever may cause instrument failure. A precondition for the use of this instrument is that an open surgical procedure should be possible as an emergency measure if the clip cannot be released from the instrument or if other unexpected circumstances arise. Always keep forceps and/or wire cutters handy to cut through the application system immediately before the kink protection if the clip cannot be released from the instrument. The use of clips to achieve haemostasis may be more difficult with hard or highly fibrotic lesions. The number of clips required for haemostasis may vary, depending on the anatomical location, histology, lesion type and the patient's condition and medical history. Do not use MRI procedures on patients in whom clips have been placed in the gastrointestinal tract. This may cause injury to patients. The use of clips in the presence of a bacterial contamination may intensify or prolong an infection. It is possible that the treatment of lesions in the oesophagus or the lesser curvature of the stomach may prove difficult. When aspirating bodily fluid, take care not to aspirate any clip that has come loose in the body cavity. This may adversely affect the aspiration function of the endoscope. Refer to the endoscope manufacturer's operating manual. If the spiral or the collet are permanently deformed and/or kinked BEFORE or DURING the use of the Clipmaster, you must STOP using the instrument. Never abruptly angulate the endoscope while the instrument is protruding from the distal end of the endoscope. This may result in injury to the patient, such as perforation, haemorrhage or mucosal injury. The use of high-frequency surgery during or after the application of clips, eg. for polyp resection, may result in injury to the patient. This is due to the risk of burns to the tissue that is in contact with the clip. Do not activate the output power until you have checked the tissue around the clip. If the spiral becomes kinked during insertion or advancement, you must stop using the Clipmaster.

**Attention with Clipmaster! If the clip sleeve is extended far enough so the next clip is visible, the system may stop functioning. If the golden sleeve tilts outside the spiral during actuation of the lever, release the clip and remove the instrument from the working channel.**


## Instructions for use

Insert the instrument into the working channel of the endoscope ( $\geq 2.8$  mm). Advance the spiral slowly and evenly and in small sections until the collet at the outlet (20 mm) protrudes from the endoscope. Avoid excessively abrupt advancement, especially in the event of resistance, to prevent any kinking of the spiral. A highly curved insertion route or an inverted endoscope may make it difficult to advance the inserted instrument, with a resulting risk of kinking at the entrance of the working channel. In this case, straighten the endoscope to allow for easier insertion. The visible part of the Clipmaster spiral between the endoscope and the handle (1) must be kept as straight as possible. The tip of the handle should point in the direction of the entrance of the working channel. Excessive bending and/or kinking of the spiral must be avoided.


Press the push button several times to place a clip. Stop as soon as the golden marking is visible.

Then pull the lever until the clip is opened.


**Attention!** Do not open the clip branches wider than  $100^\circ$ . Additional actuation of the lever will cause the clip branches to close automatically and permanently

Release the system by pressing the push button several times. Now the clip can be moved by rotating the handle in the desired direction. The clip should be as vertical as possible and placed without pressure on the tissue.

Pull the lever completely through several times to place the clip. The release of the clip can be felt and heard.

**Attention!** Pressing the clip sleeve down strongly on the tissue may cause the clip to fail.

Pull the lever completely through several times to place the clip. The release of the clip can be felt and heard.

**Attention!** If the clip is not released from the system after closing, the handle must be

separated from the remaining system before the kink protector with a gripper or a wire cutter.

If the clip CANNOT be closed from the open status, the Clipmaster can still be carefully removed from the working channel. The clip branches will fold together during removal without damaging the working channel and the endoscope.

The previous clip must be completely released before activating the next clip. The second and third clips are placed in exactly the same way as the first clip. After release of the clips the instruments can be removed from the working channel. This is done by carefully retracting the spiral.

**Attention!** Do not pull the spiral back until the clip has been completely released.

The lifetime of clips is 18 to 26 days. Subsequently applied clips drop off spontaneously by itself and are excreted naturally. Therefore applied clips don't have to be removed.

## After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

## List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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