
SPE-Series *easyGure*

SPE1-V7

404450305720

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
## Intended use

The SPE1-Y series products are accessories used with endoscopic instruments to fix guide wires during the procedure.

## Product features

## Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

## General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

Ⓢ **ATTENTION!** Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

## Contraindications

The contraindications for the SPE1-Y7 series are the same as the specific contraindications for endoscopic retrograde cholangiopancreatography.

## Possible complications

Pancreatitis, cholangitis, cholecystitis, perforation, haemorrhage, aspiration, infection, sepsis, allergic reactions to contrast agents or medicines, hypertension, respiratory depression or arrest, cardiac arrhythmia or arrest.

## Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for kinks, fractures, rough surfaces, sharp edges and protrusions. Should you discover any damages on the instruments, do NOT use them, and please inform your responsible field service contact person or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Instructions for use

### Placement of the SPE1-Y7 series on the endoscope

Push the easyCure adapter laterally onto the working channel inlet valve of the endoscope.

### Placement of the guide wire

Insert the flexible tip of the guide wire at the proximal guide wire port of the ERCP catheter or papillotome into the guide wire lumen. Push the guide wire towards the distal tip until it is flush with it. Guide the instrument into the applicable duct in accordance with the operating instructions. Release the guide wire from the proximal guide wire port to the easyCure adapter from the slotted guide wire lumen. Fix the guide wire by wrapping it around the fixing fork of the easyCure adapter.

### Placement of instruments

Introduce the placed and fixed guide wire at the distal tip of the easyCure instrument. The guide wire comes out at the distal guide wire port by bending the catheter slightly. Release the easyCure fixation while having the assistant holding it in position. Introduce the instrument through the easyCure adapter up to the distal guide wire port and fix the guide wire. Position the easyCure instrument over the fixed guide wire in short steps. Pull the instrument slowly and evenly out of the endoscope after the procedure. Make sure that the guide wire is fixed by the easyCure adapter.

### Disconnecting the guide wire

Position the guide wire and the instrument in the duct system as securely as possible under radiological control (release the guide wire fixation if necessary). Hold the instrument in the same location and pull the guide wire back under radiological control until it exits the distal guide wire port. Position the guide wire for further procedures and secure its position with the guide wire fixation. Remove the instrument.

Figure no. 1: Lateral attachment of the SPE1-Y7 to the working channel inlet valve of the endoscope:



Figure no. 2: Introduction of the instruments through the medial sealing lips of the SPE1-Y7:

Figure no. 3: Fixation of guide wires by wrapping the guide wire around the fixing fork:

## After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

### List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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