
TRP1-Serie_KOALA

TRP1-S2
TRP1-S4

404450305881

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Intended use

The TRP series products are used to retrieve endoscopically resected polyps by aspiration.

Product features

Figure 1:



Figure 2: TRP-S4 additional strainer drawers, grey

Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

Ⓢ ATTENTION! Instruments marked with this symbol are intended for single use only.

The instrument can be used immediately. Before use, check the "Use by" date on the packaging, as the instruments are only allowed to be used up to this date.

All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, sterilised or reused. Reuse, reprocessing or sterilisation can alter product properties, resulting in functional failure. Reuse, reprocessing or sterilisation furthermore entails the risk of contamination, as well as the risk of cross contamination, including the transmission of infectious diseases.

Contraindications

The contraindications for the Koala series are the same as the specific contraindications for oesophagogastroduodenoscopy, colonoscopy, sigmoidoscopy, rectoscopy and endoscopic retrograde cholangiopancreatography.

Possible complications

No precautions required.

Precautions

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Instructions for use

After removing the polyp trap from the packaging, check it for kinks and breaks. Make sure the inspection window with the suction tubing is engaged in the strainer drawer 1 position. Connect the top end of the suction tubing of the polyp trap to the suction port on the endoscope. Connect the tubing of the aspiration unit to the port on the polyp trap. Aspirate the resected polyp via the working channel of the endoscope until it is enclosed in the strainer drawer. Turn the body of the polyp trap towards strainer drawer 2 until the inspection window with the suction tubing is engaged. Remove the strainer drawer and prepare the polyp for pathological examination. Then place the strainer drawer back in the special slot on the polyp trap.

The white pre-installed strainer drawers can be swapped for the enclosed grey screen drawers (see figure 2) and used for removing further adenoids when the TRP-S4 product is used.

Evacuate other polyps via the strainer drawers, if necessary alternately. Having removed all the polyps disconnect the polyp trap from the suction port on the endoscope and from the aspiration unit.

Attention! If the inspection window with the suction tubing is in the flow-through position (viewing window outside a strainer drawer), biomaterial is conveyed directly into the aspiration unit. Remember that the strainer drawers have to be engaged in the body of the polyp trap during the surgical intervention in order to make it possible for a necessary vacuum to build up.

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

List of all icons used on FUJIFILM medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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