WIR-Series

WIR1-E1

WIR1-F1

WIR1-F5

WIR1-F7

WIR1-H9

WIR1-J2

WIR1-J6

WIR1-S2

WIR1-S6

WIR1-P1 - POL3-Series

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
Instructions for use

Intended use

WIR-series instruments are used for exploring the biliary system and guiding instruments through the barely accessible passages of the upper and lower gastrointestinal tract.

Product features

Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

ATTENTION! Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide. Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All FUJIFILM medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

FUJIFILM medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

Indications

The indications correspond to the general indications for endoscopic retrograde cholangiopancreatography, cholangioscopy and duodenoscopy.

Contraindications

The contraindications for the WIR series are the same as the specific contraindications for oesophagogastroduodenoscopy, colonoscopy, sigmoidoscopy, rectoscopy, enteroscopy and endoscopic retrograde cholangiopancreatography.

Possible complications

Pancreatitis, cholangitis, cholecystitis, perforation, haemorrhage, aspiration, infection, sepsis, allergic reactions to contrast agents or medicines, hypertension, respiratory depression or arrest, cardiac arrhythmia or arrest, foreign body in patient, inflammation, air embolism.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Do not use any alcohol, antiseptic solutions or contrast agents and solvents containing alcohol. These can damage the hydrophilic coating. Do not dry the guide wire with a cloth. Do not pull at the distal end to remove the guide wire from the dispenser. Do not touch the distal tip with your fingers.

Attention! Before the use of HF current, all guide wires except for the WIR1-S2/S6/J2/J6/H9 series must be retracted into the catheter, otherwise unintended current flow can occur under certain circumstances. Do not use any rigid and angular instruments because the coating of the guide wires could be damaged. Use disposable papillotomes exclusively to prevent damage to the lumen.

Instructions for use

After removal from the packaging, check the guide wire for kinks or breaks. Completely moisten the tip of the guide wire with sterile saline solution or sterile distilled water. Rinse the lumen of the instrument with sterile saline solution or sterile distilled water. Push the guide wire with the flexible tip forwards into the instrument. Observe the instructions for use of the instrument in this respect.

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

List of all icons used on FUJIFILM medwork products

| | | | |
|--|---------------------------------|--|------------------------------------|
| | Date of manufacture | | Use by |
| | Follow the instructions for use | | Caution near pacemakers |
| | Applied part Type BF | | Do not use if packaging is damaged |
| | Batch code | | Not suitable for lithotripsy |
| | Item number | | Contains natural rubber latex |
| | Packaging unit | | Gastroscopy |
| | Do not reuse | | Colonoscopy |
| | Sterilised using ethylene oxide | | Enteroscopy |
| | Non-sterile | | ERCP |
| | Manufacturer | | Medical device |
| | Sterile barrier system | | |

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